

One Driver for Everything, from Orthodontic Screws to Prosthetic Screws!



e-Driver

Prevent orthodontic screw fractures, secure precise insertion path, minimize screw loosening and increase the convenience of abutment fastening

All Prosthesis-related Screw Fastenings in One!

- Minimize fracture/loosening with precise fastening torque
- Easy access to posterior area by contra-angle structure
- Shorten chair time due to Rapid fastening



From Orthodontic Screw Insertion to Removal!

- 15~35Nm torque allows insertion and removal, even at hard bone
- Torque setting function minimizes fracture/detachment of orthodontic screws



e-Driver Features

1. Strong and Precise Torque

- Adjustable within the range of 5Ncm~35Ncm (5Ncm unit)
- Apply torque precisely with auto stop function
- Check torque change in real time on the LCD screen
- Precise calibration by TCS (Torque Calibration System)

2. Adjustable Speed

- Adjustable between 15RPM~60RPM (15RPM unit)
- Procedure is two times faster than manual operation
- One-touch reverse mode

3. German Technology

- Precise torque and speed implementation using high-quality German FAULHABER motor

Torque Setting Guide

Prosthesis-related Screw	Torque (Ncm)	
	Mini	Reg
Cover Screw, Healing, Imp Coping	5	5
One piece ABT (Rigid, Solid ABT)	30	30
Two piece ABT (Transfer, Angled ABT)	20	30
Temporary Abutment	20	20
Cylinder Screw	20	20

Orthodontic Screw	Torque (Ncm)	Speed (RPM)
Ø 1.4mm	15	
Ø 1.6mm	20	20 ~ 30
Ø 1.8mm	30	