ESR Full KIT

User Manual

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1. Description

ESR (Easy Screw Removal) KIT is a kit with functions for abutment & abutment screw fracture removal, slipped connection part removal, and implant’s internal screw thread restoration. The ESR KIT consists of abutment & abutment screw fracture removal device, slipped connection part removal device, and implant’s internal screw thread-restoring device.

2. Indication for use

A device intended for use in dental implant operation.

3. Instruction for Use

- Fractured screw removal guide

< Product guide by fixture type >					
Fixture type			Tool information		
Company	System	Fixture Size	Guide	SR Drill	SR Tip
Osstem Implant	TS	Mini	OGTML	OSRD08L	ORT16L
		Regular	OGTRL	OSRD10L	ORTL
	US	Mini	OGUML	OSRD08L	ORT16L
		Regular	OGURL	OSRD10L	ORTL
		Wide	OGUWL		
	SS	Regular	OGSRL	OSRD10L	ORTL
		Wide			
Nobel Biocare	Active	NP(Ø3.5)	OGNA01L	OSRD08L	ORT16L
		RP(Ø4.3)	OGNA02L	OSRD10L	ORTL
		RP(Ø5.0)			
	Replace	RP(Ø4.3)	OGNR02L	OSRD10L	ORTL
		WP(Ø5.0)	OGNR03L		
		Ø6.0	OGNR04L		
Straumann	Bone Level Implant	NC(Ø3.3)	OGSB01L	OSRD08L	ORT16L
		RC(Ø4.1)	OGSB02L	OSRD08L	ORT16L
		RC(Ø4.8)			

3i	Full OSSEOTITE Tapered Certain	Ø3.25	OGIF01L	OSRD08L	ORT16L
		Ø4.0	OGIF02L	OSRD08L	ORT16L
		Ø5.0			
		Ø6.0			
Astra	Osseo Speed TX	SMALL(Ø3.5)	OGAO01L	OSRD08L	ORT16L
		SMALL(Ø4.0)			
		LARGE(Ø4.5)	OGAO02L	OSRD10L	ORTL
		LARGE(Ø5.0)			
Zimmer	Tapered	Ø3.7	OGZB01L	OSRD09L	ORT18L
		Ø4.1			
		Ø4.7			
		Ø6.0	OGZB02L	OSRD09L	ORT18L
Biohorizons	Internal	Ø3.8	OGZB01L	OSRD09L	ORT18L
		Ø4.6			
		Ø5.8	OGZB02L	OSRD09L	ORT18L

- 1) Prepare the Guide and Screw Removal Drill (SR Drill) suitable for the fixture.
 - The Guide and SR Drill can be short (option) or long depending on the intermaxillary distance of the patient.
- 2) After connecting the Guide to the fixture, fix the Guide using an ESR Handle.
 - If the window of the Guide is connected in the direction where the operator has clear view, the operational processes are confirmed well, and irrigation to the internal space is possible.
 - In case of the external type, the Guide is more accurately fixed by using heavy impression material.
- 3) The operation engine is set to 1,500rpm and reverse rotation mode; connect the SR Drill to the hand piece.
- 4) Drilling the fractured screw by reverse-rotating the SR Drill together with irrigation.
- 5) Drill at 5~10N (less than about 1kg) while gently pumping until the red band of the handle part is gone.
 - Sometimes, a fractured screw can be hung and removed during the drilling process.
- 6) Insert the Screw Removal Tip into the Guide and remove the fractured screw by reverse-rotating it slowly.
- 7) When the internal screw thread of the fixture is damaged, restore the screw thread using Re-Tap.
- 8) Perform sufficient irrigation to remove foreign objects such as metals left inside the fixture.

• Fractured abutment removal guide

- 1) In case of 1-Piece Abutment fracture, a hole is formed using dental burr and connected by reverse rotation using the Abutment Removal Tip; the abutment is then removed by shaking using forceps.
- 2) In case of 2-Piece Abutment fracture, it is connected to the internal hole of the abutment by reverse rotation using the Abutment Removal Tip; the abutment is then removed by shaking using forceps.

• Abutment screw hex slip removal guide

- 1) In case of slip in the 1.2 Screw hex part, the hex part is slightly grinded using a small burr, and then connected using the Abutment Removal Tip for the mini option and reverse-rotated for the removal.
- 2) In case of slip in the cover screw hex connection part, it can be removed with a slot driver after making a straight “—” groove using a dental Ø0.8 burr.

• Smudging down of the internal screw thread guide

- 1) If it'll have to screw fastening fixture, There is a possibility that the thread damaged.
- 2) Prepare a tool that matches the internal information of the fixture.
- 3) Connect Re-Tap to the fixture, install a torque wrench on Re-Tap, and apply slow tapping with 30 Ncm in the normal rotating direction.
- 4) When the neck of the torque wrench is bent, remove the re-tap by reverse rotation.
- 5) Repeat normal rotation and reverse rotation until tapping is no longer possible and remove the re-tap.

4. Precautions

- Be sure to familiarize yourself with the instructions before operation.
- This product is a medical device and should be used after sterilization in the autoclave before surgery. (Temperature: 132°C, time: 15 min)
 - 1) It is recommended that the Screw Removal Drill (SR Drill) be used for not more than 1 times and with sufficient irrigation. (Do not reuse)
 - 2) The SR Drill should be connected to the Guide to prevent damage to the screw thread inside the fixture and the blade part of the drill due to shaking.
 - 3) The recommended speed of the SR Drill is 1,200~1,500rpm. Do not use at a speed of more than 2,000rpm with Low- or High-speed hand piece.
 - 4) Metal powder can be produced in the course of removing a fractured screw;

irrigate sufficiently after the operation, and then remove foreign objects using suction.

5. Maintenance and care

- Separate all used tools immediately after surgery, after cleaning and drying. Store at room temperature.
- Do not leave the instruments in a place vulnerable to contamination.
- Do not use hydrogen peroxide as a disinfectant or cleaner (it may cause damage or discoloration of the TiN Coating, laser marking and color coding).
- This product comes with a 1-year warranty from the date of purchase.

6. Side effect

- Post-operation side-effects may include local abscess, fever, temporary pain edema, hematoma, and bleeding. Qualitative and quantitative deficiency in remaining bone, allergic response, poor oral hygiene, patient in compliance, implant mobility, local tissue alteration, and improper implant positioning or arrangement may cause above complications.

7. Device information

- For device information such as part name, part number and specifications, manufacturing date, lot number, and expiration date refer to the label on the package.

* **Note** : For additional details, refer to catalogs and operation manuals. For general information on implants, refer to the catalog, prosthetic manual, surgical manual and other related data.

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DEUTSCHE OSSTEM GmbH