


OneGuide **KIT**




Keyless
open guided
surgery kit


OSSTEM[®]
IMPLANT

OneGuide^{KIT}

Keyless, precise, open guided surgery kit, designed to allow for more accurate, efficient and safe surgery with minimal error and drill heating.




- **Keyless**, open guided system.
- **Compatible** with most in-surgery software systems.
- **Tapered design** drilling sequence.
- **Easier implant placement in the posterior region.**
The keyless design allows easier drilling sequence and implant placement where patients have limited opening, allowing the drill to be inserted through the side window.
- **Reduced movement** as the drills included in the kit minimise movement during the drilling process.
- **Reduced heat at drill site.**
A side opening in the OneGuide drills facilitates a water supply, cooling the contact area to prevent heating at the drill site.
- **Shorter drilling steps.** Drilling steps are shortened dramatically by adopting 122 concepts, using different techniques according to bone density.



Precise surgery with unshaken drilling


Initial drilling:

Smart design which allows the parallel part of the drill to be in the stent for up to 3mm prior contact with the bone.



Backward drilling:

Allows precise drilling as the drill is always in contact with the guide area, due to tapered drill.




Easy access to the back of the mouth

A lateral window can be designed in the stent which facilitate drilling in the posterior region even with limited mouth opening.

Intermaxillary space:

OneGuide - 36mm,


Conventional guide - 51mm



Reduce bone heating

Enables fast drilling without heating, due to the OneGuide's side opening, facilitating douche, whilst still allowing the outstanding cutting force of the drill. (Recommended rpm: 800~1200rpm)

Outstanding cutting force drill




Side open facilitating water supply



Shortening drilling steps dramatically by adopting 122 Concepts

2 ~ 4 Drilling to implant placement protocol depending on bone density.



Bone Density	Initial	F3.5	F4.0	F4.5	Fixture
Soft	▲	▲			
Normal	■	■	■		
Hard	●	●	●	●	

Implant placement

The Digital Workflow




Data Preparation



Treatment Planning

OneGuide Guided Surgery



Digital Scanning



Prosthetic Design



Prosthetic Fabrication

OneGuide^{KIT}

This keyless guided surgery kit is designed to be simple and safe. Featuring several precise drilling tools and guide development from Hiossen® Implant, this kit allows for more accurate, efficient and safe surgery with minimal error and drill heating.



Initial Drill

- Implant location decision after Tissue punch.
- Securing the guide depth of backward drill
- Composed of 3 kinds of tools (under F4.5/only for F5.0/ only for F3.5 soft bone).




Flattening Drill

- Using for the ridge with narrow or uneven.
- Enabling Stable cutting without splattering with many of cutting edges.
- Composed of 2 kinds of tools (under F4.5/only for F5.0).



Tissue Punch

Diameter
Ø3.5, Ø4.0, Ø4.5
& Ø5.0.
Only for Tissue cutting.



NoMount Driver

F3.5: Ø5.0
F4.0/4.5: Ø5.0
F5.0: Ø5.7.
Only for OneGuide.



Fixture Driver

F3.5 F4.0/4.5
& F5.0.
Only for OneGuide.



F4.5 (Ø5.0 template hole)



Ø2.2 F3.5 F4.0 F4.5

F5.0 (Ø5.7 template hole)



F3.5 F4.5 F5.0 F5.5

OneGuide Drill

- The taper drill optimized for TSIII/IV fixture (possible to implant diameter Ø3.5~5.0, length 7~13mm fixture).
- Stable drilling with multistage structure.
- Composed of 2 kinds of tools (under F4.5, only for F5.0).

Drilling Sequence

Fixture Diameter	Bone Density	Initial	Ø2.2	F3.5	F4.0	F4.5	F5.0	F5.5	Fixture
F3.5	Soft	▲	▲						
	Normal	■		■					
	Hard	●		●	●				
F4.0	Soft	▲		▲					
	Normal	■		■	■				
F4.5	Soft	▲		▲	▲				
	Normal	■		■	■	■			
	Hard	●		●	●	●	●	●	
F5.0	Soft	▲		▲	▲	▲			
	Normal	■		■	■	■	■		
	Hard	●		●	●	●	●	●	

Implant placement